

Indications for use xT-DNA

The Tempus xT-DNA is a next generation sequencing-based in vitro diagnostic device intended for use in the detection of single nucleotide variants (SNVs) and insertion and deletion alterations (INDELs) in 648 genes using DNA isolated from formalin-fixed paraffin embedded (FFPE) tumor tissue specimens from previously diagnosed cancer patients. The xT-DNA assay is a single-site assay performed at Tempus Labs, Inc.

Limitations

The Tempus xT-DNA assay is intended to provide tumor gene alteration profiling to be used by qualified health care professionals in accordance with professional guidelines in oncology for patients with solid malignant neoplasms, however such additional tumor profiling results are not prescriptive or conclusive for labeled use of any specific therapeutic product.

Contradictions

There are no known contraindications.

Warnings and Precautions

Alterations reported may include somatic (not inherited) or germline (inherited) variants; however, the test does not distinguish between germline and somatic alterations. The test does not provide information about susceptibility to cancer. Biopsy may pose a risk to the patient when archival tissue is not available for use with the assay. The patient's physician should determine whether the patient is a candidate for biopsy.

Product Description

Tempus xT-DNA is a targeted panel of 648 genes enriched for clinically relevant genes and most commonly mutated cancer driver genes with additional genes of emerging clinical significance focused on immediately actionable mutations. Tempus works directly with oncologists and pathologists to obtain tumor samples. Once a physician and patient determine that Tempus testing is the right next step, the test can be ordered by the physician.

Our kits include everything needed for specimen collection, including easy-to-follow instructions, specimen guidelines, and packaging for convenient and fast shipping to our lab. Ordering physicians receive a comprehensive report for each patient that highlights key findings, including actionable treatments, that can be immediately translated into cancer patient care.

Ordering a Test

Tempus xT-DNA assay can be ordered using the Tempus requisition found in the Tempus CE Clinical Trials Solid Tumor Collection Kit. The requisition should be fully completed, include at least two patient identifiers, and be signed by the ordering physicians or another authorized medical professional. The Tempus CE Clinical Trials Solid Tumor Collection Kit contains all necessary instructions to package and submit samples for testing.

Results/Reports

CE Marked Results:

Tempus' xT-DNA assay is CE marked for detection of somatic single nucleotide variants (SNVs) and insertion and deletion alterations (INDELs) using DNA in solid tumors. SNV and INDEL variant types include HGVS nomenclature that describes the position and alteration; nomenclature may include coding (c.), and/or protein (p.), and/or genomic (g.) descriptors.

Professional Services:

Copy number variants (CNVs) and rearrangements detected from DNA are not included in the CE mark. Copy number variants are reported as an amplification (Copy Number Gain) or a deletion (Copy Number Loss). DNA fusions can be identified on the initial clinical report as a Chromosomal Rearrangement.

Procedure

The Tempus xT-DNA oncology assay is a 648-gene hybrid capture next generation sequencing (NGS) panel designed to detect actionable oncologic targets in solid tumor samples. This includes coverage for single nucleotide variants (SNVs) and insertion-deletion (INDEL) events using DNA extracted from formalin-fixed, paraffin embedded (FFPE) tumor samples. xT-DNA library preparation is performed using IDT unique dual-index adapters, followed by hybrid capture with custom-designed IDT xGen Lockdown probes. The final library is sequenced on an Illumina NovaSeq 6000 instrument. The clinical report provides doctors with diagnostic, prognostic, therapeutic evidence and clinical trial options based on the molecular results of the test.

xT-DNA Performance Characteristics

Validation Summary

The Tempus xT-DNA next generation sequencing assay is designed to detect actionable oncologic targets by sequencing tumor samples. The xT-DNA assay covers 648 genes spanning ~3.6mB of genomic space. From DNA sequencing, somatic insertions and deletions (INDELs) and single nucleotide variants (SNVs) are detected. The assay requires specimens with a tumor content of 20% post macrodissection. An FFPE tumor sample is sequenced, performed at 500x depth of coverage. Performance specifications are listed on Table 1 below. These results establish high sensitivity and specificity for the xT-DNA assay.

Table 1

Variant Class	Limit of Detection	Sensitivity (%)	Specificity (%)
Single Nucleotide Variants (SNV)	5% VAF	98.2	99.95
Insertions and Deletions (INDEL)	5% VAF	91.8	99.9

TEMPUS SPECIMEN GUIDELINES

To ensure accurate molecular profiling, please follow the specimen guidelines below when submitting tumor samples.

TISSUE REQUIREMENTS

- Patient material, labeled with two identifiers
- Tumor samples should be from the most recent procedure, if adequate for testing, should be less than 6 years old
- FFPE Fixation Requirements
 - 10% formalin fixation (neutral buffered) for 6–72 hours, paraffin embedded
 - No decalcification of the samples is preferred (EDTA is accepted)
- Tumor required to be at least 20%* of the sample by ratio of tumor nuclei to benign nuclei
- Optimal tumor size = 25 mm², minimum = 5 mm²

BLOCKS

*Preferred collection method



- 1 FFPE block with greatest tumor content
- 1 H&E stained slide, optional

If you would like us to evaluate quality, please send multiple blocks with a return address.

OR

SLIDES



- 10 FFPE unstained slides for NGS 5µm sections on positively charged, unbaked slides
- 1 Terminal H&E stained slide

**Submit 10 additional slides if tissue size is <25 mm²

Tempus CE Clinical Trials Solid Tumor Collection Kit

Components



FedEx shipping envelope



Biohazard bag



Sticker seal



2 Slide Cases
(each case can hold 10 slides)



Requisition form

TEMPUS xT-DNA TECHNICAL INFORMATION

Gene List

ABCB1	BLM	CHEK2**	EPH5	FGFR2	HLA-C	JAK1	MB1	PAK1	P7CH1	SH2B3	TTEC
ABCC3	BNIP3RA**	CIC	EPH82	FGFR3	HLA-DMA	JAK2	MITF	PALB2**	P7CH2	SHH	TGFB2
ABL1	BRAF	C17A	EPOR	FGFR4	HLA-DMB	JAK3	MK167	PALLD	PTEN**	SLC2A3	TGFB2
ABL2	BRC1**	CKS1B	ERBB2 (HER2)	FH**	HLA-DGA	JUN	MLH1**	PAK3	PTPN11	SLC4A2	TIGIT
ABRXAS1	BRC1A2**	CREBBP	ERBB3	FHIT	HLA-DGB	KAT5A	MLH3	PAK5	PTPN13	SLC9A3R1	TMEM107
ACTA2	BRD4	CRKL	ERBB4	FLCN**	HLA-DRI	KDM5A	MLL3	PAK7	PTPN22	SLIT2	TMEM109
ACVR1 (ALK)	BRP1**	CRLF2	ERCC1	FLT1	HLA-DPB1	KDM5C	MN1	PAK8	PTPRD	SLX4	TM6RSS2
ACVR1B	BTG1	CSF1R	ERCC2	FLT3	HLA-DPB2	KDM5D	MPL	PBRM1	PTPRT	SMAD2	TNF
AGO1	BTK	CSF3R	ERCC3	FLT4	HLA-DQA1	KDM5A	MRE11	PCBP1	QKI	SMAD9	TNFAIP3
AJUBA	BUB1B	CTC1	ERCC4	FNTB	HLA-DQA2	KDR	MS4I	PCDD1	RAC1	SMADH**	TNFRSF14
AKT1	CT1or195	CTCF	ERCC5	FOXO1	HLA-DQB1	KGAP1	MSH2**	PCDD1LG2	RAD21	SMARCA1	TNFRSF17
AKT2	C3orf70	CTLA4	ERCC6	FOXO2	HLA-DQB2	KSL	MSH3**	PDGFRA	RAD50	SMARCA4	TNFRSF9
AKT3	C8orf34	CTNNA6	ERG	FOXO3	HLA-DQA	KIF1B	MSH4**	PDGFRB	RAD51	SMARCB1	TOP1
ALK	CALR	CTNNA8	ERF1	FOXO3	HLA-DRE1	KIT	MTAP	PDK1	RAD58B	SMARCE1	TOP2A
AMER1	CARD11	CTRC	ESR1	FOXO4	HLA-DRE5	KLFF4	MTNFD2	PHF6	RAD59C**	SMC1A	TP53**
ARC**	CARM1	CUL1	ETS1	FOXO5	HLA-DRE6	KLHL6	MTNFR	PHGDH	RAD59D**	SMC3	TR6
ARL1AR	CASP8	CUL3	ETS2	FRS2	HLA-E	KLJN	MTOR	PHLPP1	RAD54L	SMO	TRAF1
APOB	CASR	CUL4A	ETV1	FUBP1	HLA-F	KMT2A	MTRR	PHLPP2	RAF1	SOC31	TRAF7
AR	CSFB	CUL4B	ETV4	FUS	HLA-G	KMT2B	MUTYH**	PHOX1B	RANBP2	SOX2	TRAF3
ARAF	CSL	CUX1	ETV5	GAPD	HNF1A	KMT2C	MYB	RAS4	RARA	SOX10	TRAF7
ARHGAP26	CBLB	CXCR4	ETV6**	GABRR4	HNF1B	KMT2D	MYC	PKSC2B	RASAF	SOX2	TSC2**
ARHGAP35	CBLC	CYLD	EWSR1	GALNT12	HOKA1	KIRAS	MYCL	PKSCA	RB1**	SOX9	TSC2**
ARID1A	CBR3	CYP19B1	EZH2	GATA1	HOKB15	L2HGDH	MYCN	PKCIB	RBM10	SPEN	TSFR
ARID1B	CCDC8	CYP2D6	FAM46C	GATA2**	HRA5	LAD3	MYO8B	PKCICD	RECQL4	SPINK1	TUSC3
ARID2	CCND1	CYP3A5	FANCA	GATA3	HSD17B2	LATS1	MYH11	PKCICG	RET**	SPOP	TYMS
ARID5B	CCND2	CYSLTR2	FANCB	GATA4	HSD3B1	LCK	NBN**	PKOR1	RHEB	SPRED1	UGAF1
ASPS	CCND3	DAX1	FANCC	GATA6	HSD3B2	LDLR	NCOB1	PKOR2	RHOA	SRC	UBE2T
ASPSCR1	CCNE1	DDI2	FANCD2	GEN1	HSP90AA1	LEF1	NCCR2	PIM1	RICTOR	SRSF2	UGT1A9
ASXL1	CD19	DDR2	FANCE	GLI1	HSPH1	LMNA	NF1	PLCG1	RNF11	STAG2	UGT1A9
ATC	CD22	DDX3X	FANCF	GLI2	IDH1	LMO1	NF2**	PLCG2	RIT1	STX3	VAMP5
ATM**	CD274 (PD-L1)	DCER1	FANCG	GNAO1	IDH2	LRP1B	NFE2L2	PML	RNF139	STX4	VEGFA
ATP7B	CD40	DRC2	FANCI	GNAO3	IDO1	LYN	NFKBIA	PMS1	RNF43	STX5A	VEGFB
ATR	CD70	DIG3	FANCL	GNAQ	IFT1	LZTR1	NHP2	PMS2**	RO51	STX19B	VHL**
ATRX	CD79A	DIGL2	FANCM	GNA5	IFT2	MADD2L2	NKX2-1	POLD1**	RPL5	STX16	VSR
AURKA	CD79B	DNAI1	FAS	GPC3	IFT3	MAF	NOP10	POLE**	RPS15	STK11**	WEE1
AURKB	CDCT3	DNAI2	FAT1	GPS2	IFNAR1	MAFB	NOTCH1	POLH	RPS6KB1	SUFU	WNK1
AXIN1	CDH1**	DNAH23A	FBIK11	GREM1	IFNAR2	MAG2	NOTCH2	POLQ	RPTOR	SUZ12	WNK2
AXIN2**	CDK12	DOT1L	FBN1	GRIN2A	IFNGR1	MAL1	NOTCH3	POF1	RRM1	SYK	WRN
AXL	CDK4	DPYD	FCGR2A	GRM3	IFNGR2	MARCK1	NOTCH4	POLQ2F2	RSF1	SYNE1	WTT**
B2M	CDK6	DYNC2H1	FCGR3A	GSTP1	IFNL3	MARCK2	NPM1	PPARA	RUNX1**	TAF1	XPA
BAP1	CDK8	EBF1	FDP5	H9	IKBKE	MARCK4	NQO1	PPARD	RUNX1T1	TANC1	XPC
BARD1	CDKN1A	ECT2L	FGF1	HGF3A	IKZF1	MARCK1	NRAS	PPARG	RORA	TAP1	XPC1
BCL10	CDKN1B	EGF	FGF10	HAS3	IL13RA	MARCK7	NRG1	PPM1D	SCDB	TAP2	XRCC1
BCL11B	CDKN1C	EGFR**	FGF14	HAVCR2	IL5	MARCK1	NSD1	PPP1R15A	SOHA	TARBP2	XRCC2
BCL2	CDKN2A**	ESLN1	FGF2	HDAC1	IL2RA	MAX	NSD2	PPP2R1A	SOHAF2**	TBC1D12	XRCC3
BCL2L1	CDKN2B	EPHA2	FGF23	HDAC2	IL6R	MCR	NTSC2	PPP2R2A	SOHB**	TBL1XR1	YEAT54
BCL2L11	CDKN2C	ELF3	FGF3	HDAC4	LTR	MCL1	NTHL1	PPP3C	SOHC**	TBK1	ZNF403
BCL6	CEBPW**	ELOC (T08B)	PGF4	HGF	ING1	MDM2	NTRK1	PRCC	SOHD**	TCF3	ZNF420
BCL9	CEP57	EMSY	PGF5	HIF1A	INP1A	MDM4	NTRK2	PRDM1	SEC23B	TCF7L2	ZNF217
BCLAF1	CFTR	ENG	PGF6	HIST1H1E	IRF1	MED12	NTRK3	PREX2	SEMA3C	TCL1A	ZNF471
BCOR	CHD2	EP300	PGF7	HIST1H3B	IRF2	MEF2B	NUDT15	PRKARIA	SETBP1	TER1**	ZNF420
BCORL1	CHD4	EPCAM**	PGF8	HIST1H4E	IRF4	MEN1**	NUP98	PRKDC	SETD2	TET2	ZNF760
BCR	CHD7	EPHA2	PGF9	HLA-A	ISG2	MET	OLIG2	PRKN	SF3B1	TFE3	ZNF573
BIRC3	CHK1	EPHA7	PGFR1	HLA-B	ITPKB	MDMT	PORYB	PRSS1	SGK1	TFEB	ZRSR2

* Includes promoter region

** Genes in which incidental germline findings are reported