

TO Validation

The Tempus Tumor Origin (TO) test uses information from analysis of nucleic acids by next-generation sequencing (NGS) performed as part of a separately-ordered Tempus|XT test. The TO test uses tumor mRNA expression results to predict the highest probability primary and other potentially likely cancer types from 68 possible diagnostic cancer types. The Tempus TO test can be reported only for samples with $\geq 20\%$ tumor purity.

Validation of the Tempus TO test was performed in a CLIA-certified, CAP-accredited lab. An analytical validation cohort consisting of 9,210 tumor samples of known origin and 1,708 cancers of unknown primary (CUPs) was created from the Tempus database. The validation cohort consisted of 25% of the labeled samples, which were selected via stratified random sampling, from within the Tempus database. All CUPs were used for the biomarker concordance analysis. All samples with an mRNA expression profile passing the quality control checks were eligible for the study. These samples were assigned one of 68 diagnostic labels or designated as CUPs based on the diagnosis assigned to the sample by Tempus pathologists at the time of sample accessioning and histologic review. Samples included formalin-fixed, paraffin-embedded (FFPE) slides, FFPE tissue blocks, blood, bone marrow aspirates, and fresh frozen tissue.

Labeled samples (n= 9,210) were used to characterize the performance of the TO test. The samples contained a broad range of tumor purities as determined by manual pathologist review, ranging from $< 20\%$ to 100% tumor purity. Performance metrics of the assay were stratified by type and subtype; the overall accuracy was 91%.

Samples labeled as a CUPs (n= 1,708) were not used in determining performance specifications of the classifier, but were utilized in other studies to characterize reportable range and interassay reproducibility.

The sixty eight (68) possible diagnostic types of the TO test are:

Acute lymphoblastic leukemia	Gastrointestinal stromal tumor	Pancreatic neuroendocrine tumor
Acute myeloid leukemia	Goblet cell adenocarcinoma	Peripheral nerve sheath tumor
Adenoid cystic carcinoma	Gynecological clear cell carcinoma	Prostate neuroendocrine carcinoma
Adrenal cortical carcinoma	Head and neck squamous cell carcinoma	Prostatic adenocarcinoma
Anogenital squamous cell carcinoma	Hepatocellular carcinoma	Renal chromophobe carcinoma
B cell lymphoma	High grade glioma	Renal clear cell carcinoma
Breast carcinoma	Leiomyosarcoma	Renal papillary carcinoma
Carcinosarcoma	Liposarcoma	Rhabdomyosarcoma
Cervical carcinoma	Low grade glioma	Salivary carcinoma
Cholangiocarcinoma	Lung adenocarcinoma	Schwannoma
Chondrosarcoma	Lung squamous cell carcinoma	Skin neuroendocrine carcinoma
Chronic lymphocytic leukemia	Medulloblastoma	Skin squamous and basal cell carcinoma
Chronic myeloid leukemia	Melanoma	Small bowel adenocarcinoma
Colorectal adenocarcinoma	Meningioma	Small cell lung carcinoma
Endometrial serous carcinoma	Mesothelioma	Synovial sarcoma
Endometrial stromal sarcoma	Metaplastic breast carcinoma	T cell lymphoma
Endometrioid carcinoma	Multiple myeloma	Thymic squamous cell carcinoma
Ependymoma	Neuroendocrine lung tumor	Thyroid cancers
Ewing sarcoma	Oligodendroglioma	Urothelial carcinoma
Fibrous sarcoma	Osteosarcoma	Urothelial neuroendocrine carcinoma
Gastroesophageal adenocarcinoma	Ovarian mucinous adenocarcinoma	Vascular sarcoma
Gastroesophageal squamous cell carcinoma	Ovarian serous carcinoma	Well differentiated gastrointestinal neuroendocrine tumor
Gastrointestinal neuroendocrine carcinoma	Pancreatic adenocarcinoma	