

Breast cancer intrinsic subtypes predict outcomes in primary and metastatic samples

Talal Ahmed, PhD¹, Mark Carty, PhD¹, Kaveri Nadhamuni, MS¹, Raphael Pelossof, PhD¹
¹Tempus Labs, Inc., Chicago, IL

INTRODUCTION

- We retrospectively analyzed 7,021 de-identified breast cancer patients with known hormone receptor (HR) or HER2 status and a matched RNASeq sample.
- The prognostic and predictive value of the PAM50 intrinsic subtypes, namely Luminal A, Luminal B, Her2, and Basal-like subtypes, is well-studied in primary as well as metastatic breast cancer settings.
- Prosigna has emerged as a rapid PAM50 subtype predictor based on the NanoString nCounter assay.
- However, assay reproducibility across various RNASeq or qRT-PCR platforms can be challenging, especially when applying the predictor on metastatic breast cancer tumors.
- Here, we create an in-house subtype predictor that works on RNA sequencing data.
- We evaluate real-world outcomes for our intrinsic subtype predictions across various immunohistochemical (IHC) labels and metastatic tumor sites.

METHODS


- The subtype predictor was trained on a cohort of 2,497 breast cancer patients using the PAM50 genes, profiled using Nanostring RNA nCounter assay (GSE148426)¹. The reference dataset contained samples collected from various sites including breast (n=2440), liver (n=936), lymph node (n=577), lung (n=540), and bone (n=304).
- For subtype classification, we first batch-corrected the external dataset to the Tempus RNA-seq reference dataset using SpinAdapt², then trained a Support-Vector Classifier (SVC) on the corrected data. A 10-fold CV experiment was performed on the corrected dataset to analytically validate the intrinsic subtype predictions.
- The concordance between HR/HER2 status and PAM50 prediction was analyzed, and these patients were excluded from training.
- Real-world overall survival (rwOS) was evaluated from the time of first diagnosis. The outcomes across PAM50 subtypes were assessed according to tumor collection site and HR/HER2 IHC status.

- Analyzing 7,021 patient samples finds a highly accurate PAM50 breast subtype predictor F1-score of 0.97 in a 10-fold CV experiment.
- The prognostic value of the basal subtype was significant for breast cancer patients across various sites of metastasis including lymph node, liver, lung, and bones.
- For patients in each of the HR+/HER2- and triple negative IHC groups, the PAM50 molecular subtypes provided an additional level of prognostic detail with statistical significance.

KEY TAKEAWAYS


RESULTS

Figure 1. Confusion Matrix with Prosigna labels and model predictions




- The 10-fold CV experiment on the Tempus-adapted GSE148426 dataset achieved F1 scores of 0.97, 0.86, 0.94, and 0.87 on Basal, HER2-like, Luminal A, and Luminal B PAM50 subtypes, respectively.

Figure 2. Confusion Matrix with IHC labels and model predictions



- Evaluating the PAM50 prediction model on Tempus Evaluation cohort, to compare the model predictions with immunohistochemical (IHC) labels.

Figure 3. Evaluation of Survival Outcomes by Cancer Type



- Basal rwOS was significantly shorter than non-basal patients (n=5,845, p< 2e-90) across all tumor sites.
- The predicted PAM50 basal rwOS remained significantly shorter than non-basal patients even when stratified by site of metastasis:
 - breast (n=2,405; p< 1e-27)
 - lymph node (n=577, p< 1e-7)
 - liver (n=936; p<1e-25)
 - lung (n=540, p< 1e-13)
 - bone (n=304, p< 1e-3).
- Among both HR+/HER2- (n=3,653) and HR-/HER2- (n=1,664) IHC cohorts with available outcomes data, the predicted PAM50 basal subtype could further stratify each of these IHC populations with basal-subtype showing significantly worse prognosis than the non-basal subtype (p< 1e-27 and p< 1e-7, respectively).

REFERENCES

- Wallden, Brett, et al. "Development and verification of the PAM50-based Prosigna breast cancer gene signature assay." BMC medical genomics (2015).
- Ahmed, Talal, et al. "Privacy preserving validation for multiomic prediction models." Briefings in Bioinformatics (2022).

Acknowledgements: We thank Vanessa M. Nepomuceno, Ph.D. for assistance with poster preparation and review.

Correspondence: talal.ahmed@tempus.com
 This presentation is the intellectual property of the author/presenter. Contact them at talal.ahmed@tempus.com here for permission to reprint and/or distribute.