# Ultrahigh Tumor Mutational Burden (TMB) Is Associated with Improved Survival Outcomes in Patients **Treated with Immune Checkpoint Inhibitors**


Douglas B. Johnson<sup>1</sup>, Binyam Yilma<sup>2</sup>, Stamatina Fragkogianni<sup>2</sup>, Metamia Ciampricotti<sup>2</sup>, Tina O'Grady<sup>2</sup>, Ben H. Park<sup>1</sup>, Justin M. Balko<sup>1</sup> <sup>1</sup>Vanderbilt University Medical Center, Nashville TN, USA; <sup>2</sup>Tempus AI, Inc, Chicago, IL, USA


Poster/Abstract #136

# INTRODUCTION

- TMB is associated with response to ICI therapy.
- A pan-tumor cutoff of 10 mutations (mt)/MB gained FDA approval of pembrolizumab although optimal cutoffs are not clearly defined.
- Here, we assessed the genomic landscape and clinical outcomes and in pts with higher cutoffs of 40 mt/MB ("ultrahigh TMB").

# METHODS






Tempus xT (DNA-seq) testing

Samples obtained within one year

before or 15 days after starting ICI



#### Analysis:

therapy


Treatment:


- Immune cell proportions (%) were estimated from RNA using quanTIseq.
- Real-world (rw) objective response rates (rwORR) were assessed at 365-days after ICI tx start, and compared using Pearson's Chi-Squared test.
- RW overall survival (rwOS) was defined as time from ICI start to death, lost to follow-up. Median rwOS was estimated using Kaplan-Meier (KM) and hazard ratios (HR) using Cox proportional hazard models.


### SUMMARY

- Ultrahigh TMB (>=40 mt/MB) comprises a novel marker for ICI response and is associated with improved rwOS independent of other prognostic markers like tumor type and MSI status.
- Ultrahigh TMB is associated with increased transcriptomic evidence of T cell and myeloid cell infiltration.


## RESULTS








rwOS was significantly improved in pts in the ultrahigh and high TMB groups compared to the low group (HR 0.38, 0.69, p<0.001; median OS x vs. y vs. z). rwPFS was also improved in ultrahigh and high TMB groups (data not shown, HR x, y, p=x).




rwOS for patients treated with ICI with various characteristics. Ultrahigh TMB resulted in superior rwOS adjusting for tumor type, MSI status, age, gender, and stage.



Characteristic	TMB Low $N = 6,699^{1}$	<b>TMB High</b> N = 2,564 <sup>1</sup>	TMB Ultrahigh $N = 475^{1}$	p-value <sup>2</sup>
B Cells	3.75 (2.79, 5.32)	3.79 (2.74, 5.39)	3.44 (2.46, 4.73)	<0.001
Nk Cells	2.37 (1.74, 3.12)	2.26 (1.72, 2.96)	2.36 (1.85, 3.05)	<0.001
Neutrophils	7.50 (5.27, 10.04)	7.66 (5.40, 10.26)	8.03 (5.70, 10.55)	0.002
Macrophages M1	5.44 (2.83, 8.94)	5.44 (2.88, 9.85)	7.31 (4.83, 10.41)	<0.001
Macrophages M2	4.01 (2.52, 6.15)	4.01 (2.52, 6.08)	3.44 (2.09, 4.92)	<0.001
Monocytes	0.00 (0.00, 0.00)	0.00 (0.00, 0.00)	0.00 (0.00, 0.00)	0.067
T Cells CD4	0.00 (0.00, 0.00)	0.00 (0.00, 0.00)	0.00 (0.00, 0.45)	0.046
T Cells CD8	0.40 (0.00, 1.35)	0.72 (0.09, 1.74)	0.83 (0.12, 1.89)	<0.001
Dendritic Cells	0.00 (0.00, 0.00)	0.00 (0.00, 0.00)	0.00 (0.00, 0.00)	0.012
Total CD4	3.88 (2.53, 6.00)	4.18 (2.73, 6.40)	4.60 (3.06, 6.43)	<0.001

Immune cell infiltration was higher in the ultrahigh TMB group for CD8, regulatory T-cells, M1 macrophages, and neutrophils, and lower for B cells and M2 macrophages (p<0.001).



Patients with higher tumor mutational burden (TMB), particularly in the ultrahigh TMB groups, showed favorable objective responses to immune checkpoint inhibitors (ICI) at 365 days. This association was more pronounced in specific cancer types like melanoma and lung cancer.

Selection Criteria

(N=17,449)