Large-scale DNA alteration prediction from H&E images across different cancer types generalizing to external images

Kshitij Ingale, Sun Hae Hong, Bolesław L. Osinski, Mina Khoshdeli, Josh Och, Rohan P. Joshi, Kunal Nagpal, Martin C. Stumpe

Tempus AI, Inc., Chicago, IL Disclosures: All authors are employees of and hold stock in Tempus Al, Inc., a for-profit company

INTRODUCTION

alterations detected by Next Generation Molecular Sequencing (NGS) have predictive, prognostic, and diagnostic value across various cancer types.

NGS may sometimes not be ordered in patients who may benefit due to low pretest probability, tissue availability, and cost considerations.

Algorithms that use digitized H&E stained images routinely acquired in the clinical workflow to rapidly screen patients for molecular alterations could be used to prioritize patients for NGS testing.

In this study, we identified a broad range of molecular alterations that can be detected using H&E images and evaluated the generalizability of a subset of these models on cases obtained from external labs and from The Cancer Genome Atlas (TCGA) projects.

METHODS

For this multi-cancer type study, a dataset of 2928 bladder cancer (and a separate dataset of 3652 with FGFR mutations), 2949 endometrial cancer, and 13,061 colorectal cancer whole slide images along with corresponding DNA alteration labels obtained from NGS was curated.

Models were trained to predict the mutation from H&E-stained images.

Additionally, logistic regression models with clinical and operational covariate features were trained as a reference.

We consider an imaging model to have a signal if it outperforms the reference model based on effect size (ROC-AUC improvement by 5%) and statistical significance (FDR corrected p > 0.05).

A subset of actionable molecular target models was evaluated on a test set consisting of slides stained at external labs as well as a test set with slides sourced from TCGA projects.


SUMMARY

- strategy for patient prioritization.
- real-world deployment.

RESULTS

Out of 78 imaging models, 46 (59%) outperformed corresponding reference models as shown in **Figure 1**, suggesting an association of morphological features with many DNA alterations. The model performance is stable on data obtained from external labs with a mean drop of 0.01 in ROC-AUC (Figure 2a). The data from TCGA projects is obtained from a diverse set of external labs and clinical characteristics, however, the models appear to be reasonably robust to these cases as well, resulting in a mean drop of 0.04 in ROC-AUC (Figure 2b).

Figure 1. DNA alterations can be detected using H&E images



• Our work suggests that H&E image-based identification of actionable molecular mutations can be a viable

• Some of the actionable target algorithms were also generalizing to external, demonstrating potential for

• Further rigorous model development can be explored for integrating these models into clinical practice.

Figure 2. Model generalizes to externally stained slides and TCGA slides